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©Copyright owned by Qingdao ChuangMeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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The Model 1104 digital display rotarional viscometer is a new type of viscometer designed by our company. The viscometer uses embedded structure design, built-in high-precision angle sensor, stepper motor drive, speed more stable and accurate; panel keys clearly marked for easy operation, high-resolution liquid crystal display can real-time display the viscosity of the test solution and stirring speed.

This product can be widely used in petroleum, chemical industry, food, medicine, cosmetics and other industries and scientific research units. Working Principle: Put the slurry cup with mud on the pallet. The mud enters the annular space (shearing clearance) formed between the outer sleeve and the float of the viscometer. When the stepper motor drives the outer sleeve to rotate, the sleeve will react with viscous resistance in the mud. This force produces torque to the float in the sleeve through the mud in the annular space, and the float is connected to a set of calibration components. When the torque is in equilibrium with the viscous resistance, the dial assembly will stabilize at a certain scale. Because the scale value is proportional to the viscous resistance of the sleeve, the final viscosity value can be obtained by multiplying the scale value by a specific coefficient and displayed on the LCD screen for reading.

There is a linear relationship between shear rate and viscosity. So

the relationship between shear stress and shear rate is linear. In practice, many fluids may not follow Newton's law, but their rheology is close to Newton's law, so viscometers can still be used to measure, and the measured viscosity is more accurate. The calibration method of model 1104 Digital Rotating viscometer is linear calibration in Newton mode. It is not suitable for the samples which do not conform to Newton linear calculation method. In this case, the nonlinear calculation method should be used for viscosity measurement and velocity calculation.

When the standard outer sleeve, floater and torsional spring rotate in 300 rounds/minute, the measurement unit of the viscometer is centipoise (or milli pascal/second). At other rotating speeds, the read data must be conducted a series of simple multiplying conversion. The sixth part will present the calculation method of viscosity of plastic fluids (such as drilling fluid).

If choosing different rotating speeds or using different combinations of torsional spring-floater, the scope of shearing rate may have some changes. In order to expand the scope of shearing stress and thus measure all sorts of liquids, we have designed different series of torsional springs, which can carry out reciprocal conversion easily.

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
Structure diagram of digital display rotary viscometer

序号 No	名称 Name	序号 No	名称 Name	序号 No	名称 Name
1	底座护罩 Base cover	9	转筒 Rotor	17	主动同步带轮 Active synchronous pulley
2	驱动器 Driver	10	刻度盘组件 Dial assembly	18	电机后护罩 Rear shield
3	托板 Supporting plate	11	从动同步带轮 Driven synchronous pulley	19	电机 Electric machinery
4	星形把手 Star grip knob	12	测量弹簧组件 Measuring spring assembly	20	"O"型圈 "O"-ring
5	外套筒 Outer sleeve	13	电机前护罩 Front shield	21	立柱 Upright column
6	内转筒 Floater	14	显示器线路板 Display circuit board	22	船型开关 Rocker switch
7	液位线 Liquid level line	15	控制面板 Control panel	23	DC 电源插座 DC power socket
8	浆杯 Slurry cup	16	同步带 Synchronous pulley		

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Warning: Operators must be familiar with the operation process of the instrument and understand potential hazards. The instrument can be supplied from 100 V to 240 V. Be sure to keep your hands, clothes and other articles away from the rotating parts of the instrument.

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The viscometer can be heated by electric heater. Before heating, ensure that the heating body is fully contacted with the sample to be tested, and ensure that the ground is grounded.

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When performing maintenance, be sure to turn off the viscometer and unplug the power. If the liquid to be tested overflows, please wipe it clean with a soft cloth. It is strictly prohibited to clean directly with

water and pour water into the machine base to avoid damaging electrical components.

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This viscometer cannot test samples with temperatures above 93 °C.

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When using a heater to heat the test liquid, be sure to prevent the test liquid from splashing out due to heat. It is strictly prohibited to touch the heater directly with hands. The heating temperature of the heater shall not exceed 93 °C.

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
The liquid is injected into the slurry cup to make the liquid reach the scale line (350ml) on the cup body. Place the slurry cup on the tray and the liquid to be measured enters the rotary sleeve, where the level line on the sleeve can indicate the immersion depth. If the immersion depth exceeds the maximum value, it may cause damage to the bearing of the float. If other sample containers are used, the distance between the bottom of the outer sleeve and the bottom of the container should be greater than or equal to 1.27cm.

Model 1104 digital rotary viscometer has twelve rotational speeds, ranging from 1 to 600 rpm. When used, the power supply is switched on first, and then the required speed is selected by the digital key on the control panel (Fig. 2). The speed and viscosity can be real-time viewed through the LCD screen.

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Outer sleeve, float and torsion spring

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
The R1-B1-F1 combination is suitable for all types of viscometers, and other combinations may be used to calculate the shear rate of the sample to be measured. When some combinations make the reading of shear stress more erroneous, the combinations can not be selected if they do not meet the requirements.

1. 外套筒的拆除与安装 Removal and installation of outer sleeve

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According to figure III, the outer sleeve can be loaded and unloaded.

Attention should be placed to the highest position when installing.

图（3）外套筒的拆装 Removal and installation of outer sleeve

2. 浮子的拆除与安装 Removal and installation of float

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When the float is removed, the float is pulled down and rotated counterclockwise; when the float is installed, the float is pushed up and rotated counterclockwise.


3. 扭力弹簧的拆装 Removal and installation of torsion spring

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When the spring is removed, the bolts C and D can be removed after removing them with reference to Fig. IV. When the spring is installed, the bolts C and D can be tightened after the spring is positioned. Be careful not to stretch the spring to prevent spring damage from exceeding the elastic limit. Before screws C and D are tightened, press or pull the spring gently to ensure that the top of the adjusting mandrel A is higher than the top of the spring knob B. The screw D should be fixed together with the spring knob B.


Note: the E shown in the following figure is a high-precision angle sensor circuit board. Attention should be paid to protection during assembly and disassembly.

图（IV）扭力弹簧的拆装 Removal and installation of torsion spring

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The viscometer shall be regularly maintained and calibrated by applying torque to the float shaft. Through timely maintenance and calibration, problems can be found in time, so that the instrument can maintain accurate measurement accuracy. There are two ways to calibrate:

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Static load calibration method is easy to operate, when the spring is adjusted, the adjustment results can be easily confirmed. The static load calibration tool is NLJ-A type torsion spring dynamometer (Figure V). The following is the calibration process description:

1) Remove the outer sleeve, the float will rotate counterclockwise (Axis taper coordination of float and float) and push upward.

2) Take a piece of inelastic filament and fix one end of the filament on the surface of the float with adhesive paper. Then wrap the filament counter-clockwise around the surface of the float for 2-3 turns. The other end of the filament is hooked on the end after passing through a special force measuring frame or a horizontal fixing wheel.

3) Select appropriate weight according to table III.

4) Adjust torsion springs as required. Reference can be made to the C "Disassembly and Installation of Torsional Springs" section of "VI. Outer sleeves, floats and torsional springs".

There are two kinds of tolerances for F1 spring matched with this

instrument: $127 \pm 0.5 \text{ mPa} \cdot \text{s}$ for 50g and $254 \pm 0.5 \text{ mPa} \cdot \text{s}$ for 100g. 0.5mP. s error is allowed when spindle rotates. This error usually decreases when testing fluid. At least read three times and get the average value. If the spring does not change linearly, it indicates that the bending of the float shaft needs to be repaired.

扭力弹簧总成 Torsion spring assembly (R1-B1 组合)	扭力弹簧常数 Torsion spring constant K1 Dynes/cm/def	扭力弹簧系数 Torsion spring coefficient	重量 (克) Weight (g)				
			10	20	50	100	200
			粘度值 Viscosity value				
F-0.2	77.2	0.2	127.0	254.0	—	—	—
F-0.5	193.0	0.5	50.8	101.6	254.0	—	—
F-1	386.0	1	25.4	50.8	127.0	254.0	—
F-2	772.0	2	—	25.4	63.5	127.0	254.0
F-3	1158.0	3	—	—	43.0	84.7	169.4
F-4	1544.0	4	—	—	—	63.5	127.0
F-5	1930.0	5	—	—	—	50.8	101.6
F-10	3860.0	6	—	—	—	—	50.8

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Torsion spring stiffness linear test table (assumed radius 1cm)

砝码(克) Weights (g)	允许转动范围(格) Allowable rotation range	砝码(克) Weights (g)	允许转动范围(格) Allowable rotation range
5	21.55~22.21	40	172.45~177.70
10	43.11~44.42	45	194.00~199.90
15	64.66~66.63	50	215.57~222.13
20	86.23~88.85	55	237.12~244.30
25	107.78~111.06	60	258.68~266.56
30	129.34~133.28	65	280.23~288.77
35	150.89~155.49	70	301.80~310.99

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F0.2 torsion spring stiffness linear test table (radius 1.725cm)

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This procedure is to be used for calibration using only Newtonian certified calibration fluids. Calibration Fluids are available in nominal 20, 50, 100, 200, and 500 cP. All are traceable to ASTM standards and each bottle of fluid is furnished with a viscosity temperature chart certifying that batch of fluid.

1. Before immersing the outer sleeve and float into the standard solution, ensure that the instrument to be tested is clean. If necessary, remove the outer sleeve and wash the float thoroughly. Make sure the float shaft and outer sleeve are intact.
2. Add the calibration fluid to the liquid level line of the slurry cup and place the slurry cup on the tray of the instrument. Lift the tray upwards until it reaches the proper depth. Reference figure I
3. Place the thermometer in the tested sample and choose a safe position to prevent breakage.
4. Operate the instrument at 300 rpm for three minutes, This will equalize the temperature of the float, rotor and the fluid.
- 5 .Record the viscosity readings on the LCD display at 300 rpm and 600 rpm. The reading of the thermometer is accurate to 0.1°C (0.15°F) .

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Newtonian Viscosity in centipoise may be read directly from the dial when viscometer is run at 300 rpm with R1-B1-F1 combination. Other springs may be used providing the dial reading is multiplied by the "f" factor (spring constant).

To rapidly determine Newtonian viscosities in cP with FANN viscometer, use the following formula:

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1 mPa.s 等于 1 cP

Example:

Using an R2-B1 combination at a speed of 600 rpm with an F5.0 spring, and a dial deflects to 189.

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Note:

Combinations with the larger gaps are likely to give results that differ from these figures. For best accuracy, calibrate with a standard fluid having a viscosity near the range of interest and using the R-B-F combination to be used in the test

外套筒-浮子组合 Outer sleeve --float combination	外套筒-浮子因子 Outer sleeve --float factor
R1-B1	1.000
R1-B2	8.915
R1-B3	25.392
R1-B4	50.787
R2-B1	0.315
R2-B2	8.229
R2-B3	24.707
R2-B4	49.412
R3-B1	4.517
R3-B2	12.431
R3-B3	28.909
R3-B4	57.815

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$$K1 = 50 \times 1.725 \times 981 / 386 = 219.2$$

Example:

Set F1 torsional spring, spring constant is 386 Dynes/cm/° def and R1

and B1 combination. Use 50 gram weight. The formula is:

$$K1 = 50 \times 1.725 \times 981 / 386 = 219.2$$

C. 数据测试及计算 Data testing and calculation:

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Adjust the temperature indoors and make it $20 \pm 5^\circ\text{C}$. Operate it according to testing operation methods. If you measure in well site, please reduce the sampling time and record sampling place and conditions on the measurement table.

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1. Newtonian fluid Absolute viscosity:

Adjust rotary speed of the instrument to 300 r/min. When the reading on the dial is constant, the reading is absolute viscosity.

$$\eta = 300\text{r/min (reading)} \text{ mPa} \cdot \text{s}$$

2. Plastic viscosity:

1) Adjust rotary speed to 600 r/min. When reading on the dial is constant, half of reading is apparent viscosity.

2) Adjust rotary speed of the instrument to 300 r/min. The difference of 300 r/min and 600r/min is plastic viscosity.

3) Mix the drilling fluid 10 minutes at high speed. Rotate at 3 r/min

then maximum reach is initial gel strength. After 10 minutes standing and then you can record get strength.

apparent viscosity: $\eta_a = 1/2 \times 600 \text{ r/min (reading)} \quad \text{mPa} \cdot \text{s}$

Plastic viscosity: $\eta_p = 600 \text{ r/min (reading)} - 300 \text{ r/min (reading)} \quad \text{mPa} \cdot \text{s}$

yield point: $\tau_0 = 0.511 (300 \text{ r/min reading} - \eta_p) \text{ Pa}$

gel strength: $\tau_{\text{begin}} = 0.511 \times 3 \text{ r/min (reading)} \quad \text{Pa} \quad (1 \text{ min})$

$\tau_{\text{end}} = 0.511 \times 3 \text{ r/min (reading)} \quad \text{Pa} \quad (10 \text{ mins})$

3. Pseudoplastic fluid:

Its flowing characteristics is that viscosity decreases with shear stress increasing. The flowing of pseudoplastic fluid follows power function.

$$\tau = k \left(\frac{dv}{dx} \right)^n \quad \lg \tau = \lg k + n \lg \frac{dv}{dx}$$

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
Production address: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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After inspection, Qingdao Chuangmeng Instrument Co., Ltd The products produced meet the requirements of the above standards. Approved for delivery.

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co., Ltd.