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©Copyright owned by Qingdao ChuangMeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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The portable roller heating furnace produced by Chuangmeng Instrument is equipped with an imported precision controller and adopts the latest microcomputer intelligent control technology. It is directly set, digitally displayed, and can provide deviation indication. The machine runs smoothly, with low noise, high precision, wide working range, high degree of automation, simple operation, sturdy and durable, and easy maintenance.

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
1. Box body: Made of steel frame structure, aluminum silicate insulation layer, stainless steel shell, sturdy and durable.

2. Drum: Made of high-quality metal material drum and frame, it is lightweight and rotates smoothly.

3. Heating control: PID intelligent control is adopted. The test time can be set and an automatic alarm can be triggered.

4. Power system: It is output by a reduction motor and driven by a sprocket to rotate three rollers. The transmission is smooth and reliable.

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The working principle and performance are fully designed according to American API standards. Simulate the slow circulation speed of drilling fluid in the downhole circulation space and study its heating and rolling characteristics.

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1. Confirm that the voltage matches and the power supply is reliably grounded.

2. Connect the power and turn on the power switch. The current temperature inside the heating furnace is displayed above the temperature controller display window, and the set temperature is displayed below. Start the motor switch, and the motor drives the drum inside the furnace to rotate, allowing for sample mixing. Make the motor idle for minutes, and the sound and speed are normal. After the preparation work is completed, normal work can proceed.

3. Temperature setting: Turn on the heating button, and the actual temperature inside the roller furnace will be displayed above the window.

Press the "∧, ∨" keys on the thermostat to set the desired control temperature at the bottom of the window. The maximum working temperature of the roller furnace is 180 °C. When the temperature approaches the set temperature, the system enters the automatic constant temperature control state.

Warning: Setting keys "<<PF" "" "" are strictly prohibited for non factory technicians to touch.

4. Timer setting: Open the cover under the timer, press the "+" or "-" keys to set the timer time ("H" for hours, "M" for minutes, "S" for seconds), start the timer switch, and the machine will automatically stop heating when the set time is reached.

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1. When moving or maintaining instruments. Handle with care to avoid deformation of components that may affect accuracy and usability.

2. After each use, the inner cavity of the roller furnace should be cleaned promptly and placed in a dry environment.

3. Regularly lubricate the sprocket and chain.

七、故障的判定与排除 Determination and elimination of faults

故障 Fault	原因 Reason	维修方法 Maintenance Method
接通电源，电机转动，滚筒不转 Connect the power, the motor rotates, but the drum does not rotate	①链轮与轴连接的紧定螺钉松动 ②链条与链轮脱开 ① The fastening screw connecting the sprocket and shaft is loose ② Chain and sprocket disconnect	①将所有链轮与轴连接的紧定螺钉拧紧 ②将脱开的链条重新按装于链轮上 ① Tighten all the fastening screws connecting the sprockets to the shaft ② Reinstall the detached chain onto the sprocket
接通电源，仪器无电 Connect the power supply, but the instrument has no power	1、电源插头接触不良 2、保险故障 1. Poor contact of power plug 2. Insurance malfunction	①检查电源插头是否牢固 ② 更换保险 1、Check if the power plug is secure 2、Replace insurance

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Production address: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.