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©Copyright owned by Qingdao ChuangMeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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The performance of drilling fluid is not only related to its components and their interactions. And it is related to the mixing device, mixing method, and conditions. Therefore, in order to obtain a correct inspection result, it is necessary to standardize all the mixing equipment, mixing methods, and conditions.

The flip type slurry mixing machine produced by Chuangmeng Instrument is a standardized slurry mixing instrument widely used in oilfield laboratories, geology, metallurgy, as well as universities and research institutes.

Adopting a flip structure design, it is flexible and convenient to use. The spindle speed adopts stepless speed regulation, with stable speed and large starting torque. The control part adopts advanced digital control circuits to ensure the smoothness of its characteristic working curve and minimize its influence from the stirring medium.

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
Adopting a flip design for infinite speed regulation, it has a small size and is easy to transport. The shell and the entire system are made of high-quality stainless steel materials that are durable and

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Pour the sample into the bucket according to the specified ratio, stir it through the mixing bucket, and adjust the stirring speed to ensure that the density of each part of the liquid is the same while fully stirring.

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
Pour the prepared sample into the mixing bucket (as shown in Fig 1).

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Flip the motor up to the stop position, place the mixing drum in the base for rotation and fixation, then reverse the direction of rotation to restore the motor to its original position, connect the power supply, press the switch, rotate the speed control knob, and press the switch again to stop the mixing motor from working. The mixing motor can adjust the mixing speed as needed at any time during operation.

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After the mixing is completed, turn off the power after the mixing shaft has completely stopped, flip the motor to the stop position, remove the mixing drum, and reset the flipping motor. (Fig 3)

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After the mixing is completed. The power should be unplugged and the machine should be wiped in a timely manner.

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Regularly check if electrical components are functioning properly.

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操作步骤如下 The operation steps are as follows:

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The experimenter is familiar with the entire operation process and the possible situations that may occur during the operation.

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Open the base of the fixed mixing drum and check if the electrical components are functioning properly.

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Open the motor cover and check if the motor and wiring are normal.

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Check if the connecting nut in Fig 1 is loose.

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When repairing and moving instruments, the power should be cut off in order to troubleshoot and move them.

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The mixing shaft must not collide to avoid bending and deformation, which may damage the instrument.

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After each use, the instrument should be promptly wiped clean and placed in a dry environment.

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Production address: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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结论: Conclusion:

经检验, 青岛创梦仪器有限公司生产的产品符合上述标准的要求。准予出厂。
After inspection, Qingdao Chuangmeng Instrument Co., Ltd The products produced meet the requirements of the above standards. Approved for delivery.

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.