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©Copyright owned by Qingdao ChuangMeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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The capillary absorption time measuring instrument produced by Chuangmeng Instrument can be used in multiple fields, such as sewage treatment (rapid testing of sludge filtration characteristics and status), testing and research of oilfield chemical treatment agents (evaluating the ability of mud to inhibit dispersion, anti swelling ability of drilling fluid, inhibition of hydration dispersion ability of drilling fluid, optimizing the formula of the working fluid entering the well, such as the type, dosage and ratio of treatment agents, etc., and can also be used as a new method for analyzing and evaluating shale classification).

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It includes a timer, stainless steel cylinder, test probe, etc. The stainless steel cylinder has a diameter of 1 centimeter at one end and 1.8 centimeters at the other end.

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The capillary absorption time analyzer measures the time required for various test solutions and shale powder slurries to permeate a certain distance through specially designed filter paper. This value is called the CST value. Its size is related to factors such as the properties of the liquid and the dispersibility of colloids, and can be used to determine the degree of colloidal dispersion of shale in water. The smaller the CST value, the better the inhibitory effect, and its minimum value indicates that:

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The minimum shale hydration effect;

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Collect rock samples, preferably with a diameter greater than 6.4mm.

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Clean the rock debris with fresh water.

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Strongly stir the rock debris that has been cleaned with fresh water using a 3% (mass fraction) H_2O_2 solution.

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Dry the cleaned rock fragments in a constant temperature drying oven at $(105 \pm 3)^\circ C$ for 2 hours.

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Grind the dried sample with a grinder.

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After passing the sample through a 100 mesh standard sieve, put it into a wide mouthed bottle for later use.

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Quantitatively weigh 7.5g of shale sample passing through a 100 mesh sieve, pour it into a stainless steel cup, and add distilled water to 50ml.

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Place the cup containing the sample on the corrugated mixer and stir at 3 speeds for 20 seconds.

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Take out 3ml of slurry using a syringe without a needle (5ml) and press it into a CST cylindrical slurry testing container.

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Continue stirring the remaining slurry at three different speeds for 60 seconds and 120 seconds respectively, and measure its CST value (i.e. repeat steps 3 and 4 in each section).

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Use 20s, 60s, and 120s as X values, corresponding CST values as Y values, and substitute them into linear regression.

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Calculate the slope m , intercept b , and correlation coefficient R in $Y=mX+b$, and plot the CST-T curve. T is the stirring time. In the regression equation $Y=mX+b$, Y represents the CST value, m represents the slope of the line, and the CST value varies with shear time, which can be used to characterize the speed of hydration dispersion. B is the intercept on the CST axis, which is derived from the extension of the CST-T curve. 一般认为, b 越大, 瞬时破裂下来的胶体颗粒越多; m 越大, 水化分散的速度越快, 反之亦然。

It is generally believed that the larger b , the more colloidal particles will break down instantly; The larger the m , the faster the hydration and dispersion speed, and vice versa.

(3) 泥浆的分析、评价和分类 Analysis, evaluation, and classification of mud

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The CST instrument can be used as an instrument for analyzing, evaluating, and proposing new methods for shale classification, as shown in the figure. It is used to determine the dispersibility of shale. The CST device consists of a filtering funnel (a stainless steel cylinder with a diameter of about 2.54cm and a height of 5-6cm), standard porosity filter paper, a timer, and electrodes connected to it. The distance

between the electrodes and the edge of the funnel is 0.5cm, 1.0cm, and 1.5cm, respectively.

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The shale dispersion test procedure for CST device is as follows: 15% of 100 mesh magma solution is sheared at a constant speed for different times, the funnel is placed on standard filter paper, and the filter paper covers the charged electrode of the test plate. Take 5ml of stirred slurry and pour it into a funnel. Measure the time required for the slurry to flow on filter paper for a distance of 0.5cm. The same experiment should be conducted at least three times, with an error not exceeding 3% ~ 5%.

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In order to evaluate the effect of different electrolytes on shale, when preparing CST test samples, it is required to rinse shale rock fragments with distilled water until there are no chloride ions present in the water, and then dry the rock fragments.

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The CST dispersion test curves of three types of shale are shown in Figure 2. The CST value in the figure is linearly related to the shear time, therefore, the shale dispersion equation can be expressed as:

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Where: y - the time required for the slurry to penetrate 0.5cm, for example, s; $M1$ - Slope, representing the dispersion rate of shale in solution; X - shear time, s; B - intercept, represents the instantaneous fine dispersion of colloidal particles (initial dispersion).

(4) 页岩的分析、评价和分类 Analysis, evaluation, and classification of shale

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The CST (Capillary Filtration Timer) test is a fast and reliable method for demonstrating the filtration performance and state of sludge, and it has been widely used abroad since 1970. The filtration performance of sludge affects the output of almost all types of dewatering equipment, including drying beds, belt presses, vacuum filters, filter presses, and centrifugal separators.

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The capillary suction pressure generated by standard filter paper is used to "absorb" water from the sludge, and the rate of water infiltration into the filter paper depends on the state of the sludge

and the filtration properties of the sludge cake formed on the filter paper. CST is obtained from two electrodes placed at standard intervals on a funnel. The time obtained from the beginning to passing between these two electrodes is called CST.

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The force generated by capillary water absorption in the funnel is much greater than the peak of fluid static pressure, so the test does not depend on the amount of sludge. As long as there is sludge, CST can be produced.

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Clean and dry all components for later use, and place the instrument in a dry environment.

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When moving or maintaining instruments. Handle with care to avoid deformation of components that may affect accuracy and usability.

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Production address: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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After inspection, Qingdao Chuangmeng Instrument Co., Ltd The products produced meet the requirements of the above standards. Approved for delivery.

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.